



# HUGOPHARM

*Systems Engineered with Mind and Spirit*

**Hugopharm Technologies Pvt. Ltd.**

Regd Office: 8, Jogani Industrial  
Estate, 541 Senapati Bapat  
Marg, Dadar (W), Mumbai 400028

Contact : Tel: 91-22-24226882 / 91-22-  
24222815 ; sales@sbpanchal.com

**Our Ref No: HUGO/2301/2025-2**

**23-01-25**

To,

|                         |                                                           |
|-------------------------|-----------------------------------------------------------|
| <b>Contact</b>          | Dr. Tarak Mehta                                           |
| <b>Business Address</b> | <b>Aculife Healthcare pvt. Ltd.</b><br>Ahmedabad, Gujarat |
| <b>Business Phone</b>   | +91 9986553811                                            |
| <b>E-mail</b>           | tarakmehta@aculife.co.in                                  |

*We are glad to quote for the following R&D formulation line of equipment manufactured by us.*

- 1. UNIPEL100- PILOT AUTOMATIC PELLETIZATION SYSTEM**

Sincerely,

**For Hugopharm Technologies Pvt. Ltd.**  
Hiren Panchal

Tel: 022-24226882 / 022-24222815

## 1. UNIPEL-100: PILOT AUTOMATIC PELLETIZATION SYSTEM

| <i>Sn.</i>                     | <i>Item</i>                     | <i>Description</i>                                                                                                                                                                                                                                                                                                                                                                |
|--------------------------------|---------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>USSE100: SCREW EXTRUDER</b> |                                 |                                                                                                                                                                                                                                                                                                                                                                                   |
| 1.                             | <b>System Description</b>       | <p><b>UNIEXTRUDE100</b>– Screw Extruder for processing extrudates 0.5mm – 3mm size.</p> <p>This method of extrusion has single screws transferring the material to the pressing cams, which in turn press the material on the required size of perforated mesh. The extruder shall have a conical discharge assembly.</p>                                                         |
| 2.                             | <b>Capacity</b>                 | <p>UNIEXTRUDE - 100 Maximum processing capacity: 15 Kg/hr.</p> <p>Note : Output depends upon the Bulk density and material specific properties. Maximum output mentioned is for 1mm Extrudate with MCC Lactose Formulation</p>                                                                                                                                                    |
| 3.                             | <b>Material Of Construction</b> | <p>All the contact parts shall be made in S.S. 316 quality material. All non contact parts shall be made from S.S. 304 quality material except bought out materials like motors and gear box.</p>                                                                                                                                                                                 |
| 4.                             | <b>Extruder</b>                 | <p>Extruder assembly consists of a conveying screw &amp; pressing flute. The assembly is closed by quick remove tri-clover clamps and hinge design for ease in operation. Standard cone mesh with 1mm die shall be provided with the machine. The barrel shall be jacketed for recirculation of chilled water for cooling.</p> <p>Note: Chilled Water supply in client scope.</p> |
| 5.                             | <b>Feeding Assembly</b>         | <p>25liter Capacity feed hopper mounted above the extrusion assembly; the function of this hopper is to feed the transfer screws. The feeding hopper shall be made from S.S 316 quality material It shall have a delumper installed in it along with speed control of the feed screw to allow dosing of the wet granules to the extruder on volumetric basis.</p>                 |

|                                 |                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|---------------------------------|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.                              | <b>Gear Box And Motor Assembly</b> | <p>A suitable gear box consisting of a pair of Spur Gear &amp; Worm gears shall be used to transmit power.</p> <p>A variable speed drive of ABB OR EQUIVALENT MAKE shall be used to regulate the speed of the screw from 20-100 R.P.M</p> <p>3 HP TEFC, 1440 RPM, IP-55 flange mounted, voltage and Hz as per customers requirement.</p>                                                                                                                                                                                                                                                         |
| 7.                              | <b>Hold Up Vessel</b>              | <p>Hold up vessel with discharge valve to shall be located between the extruder and spheronizer for feeding the spheronizer. It shall be made from S.S. 316 quality steel.</p>                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>UNISPHERE500-SPHERONIZER</b> |                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| 8.                              | <b>System Description</b>          | <p>UNISPHERE500-Spheronizer generates spherical shape pellets which come out from an extruder. It and process downstream of any extruder.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| 9.                              | <b>Capacity</b>                    | <p>USP-500 Maximum processing capacity : 1kg to 3kg/batch.</p> <p>Note : Output depends upon the Bulk density and material specific properties.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 10.                             | <b>Material of Construction</b>    | <p>All the contact parts shall be made in S.S. 316 quality material. All non contact parts shall be made from S.S. 304 quality material except bought out materials like motors and gear box</p>                                                                                                                                                                                                                                                                                                                                                                                                 |
| 11.                             | <b>DRUM ASSEMBLY</b>               | <p>A 500mm drum will be of cylindrical design to enable discharge on any side of the drum. It shall be fixed with dynamically balanced chequered plate having 3.25mm pitch. 5 HP TEFC, 1500 RPM ELECTRIC MOTOR shall be provided to rotate the friction plate.</p> <p>It shall be easily removable for cleaning with screw arrangement.</p> <p>Drum cover. Made in Acrylic. It shall have a small opening for feeding material during operation.</p> <p><b>The drum shall be jacketed for recirculation of chilled water for cooling.</b></p> <p>Note: Chilled Water supply in client scope.</p> |
| 12.                             | <b>DISCHARGE ASSEMBLY</b>          | <p>Sealed leak tight fitting, Pneumatic cylinder fixed with tubings for discharge. Hopper made in Stainless Steel 304 fixed with acrylic cover for discharge.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                |
| 13.                             | <b>Controls</b>                    | <p><b>For Common PLC System between Extruder and Spheronizer with Mitsubishi PLC and 7" Touchscreen.</b></p> <p><b>The PLC system shall control the extruder feeder, Extruder, Hold Vessel, Spheronizer. The client shall be able to control the spheronization process with better control through the PLC System as the input to each process shall be consistent.</b></p>                                                                                                                                                                                                                     |

## PRICE: UNIPEL 100: PILOT PELLETIZATION PLANT

| Sn. | Item                                                                                                                      | Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-----|---------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.  | Rs. 40,00,000/-                                                                                                           | For Automatic Pelletization Plant as described above comprising of the following 1. Screw Extruder: Model USSE110<br>2. Extuder Cone Mesh and Spheronizer Plate suitable for 18-25 Mesh pellets<br>3. Feeder for Extruder<br>4. Stainless Steel Platform for the System<br>5. Hold UP Vessel<br>6. Mitsubishi PLC with Touchscreen System.                                                                                                                                                                                                                                                 |
| 7.  | Rs. 2,00,000/-                                                                                                            | For Extruder Mesh with 2mm Size                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 1.  | Rs. 2,50,000/-                                                                                                            | For Extruder Mesh 0.8mm Size                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| 2.  | Rs. 3,00,000/-                                                                                                            | For Extruder Mesh with 0.5mm Size.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| 3.  | Rs. 2,50,000/-                                                                                                            | For Spheronizer Plate for pelletization of spheres for 2mm Pitch                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 4.  | Rs. 3,00,000/-                                                                                                            | For Spheronizer Plate for pelletization of spheres for 1mm Pitch                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| 8.  | Rs. 1,25,000/- for maximum of 2 persons for maximum 5 days.<br>Each additional days shall be charged at Rs. 5,500/- extra | <b>For Installation &amp; Commissioning.</b><br><b>CLIENT'S SCOPE OF SUPPLY:</b><br>1. Level floor.<br>2. Foundation and necessary masonry work if required.<br>3. Shifting of equipment at site.<br>4. Electrical wiring upto the equipment as per our instructions.<br>5. Steam connection along with necessary steam fittings such as steam trap, strainer etc. and steam controlling device.<br>6. Necessary unskilled help at your cost.<br><br><b>TO and Fro Travel Charges from Mumbai to Your Site, local lodging, boarding and conveyance shall be in Clients Scope of supply</b> |
| 9.  | Rs. 55,000/-                                                                                                              | For Set of Validation Documents comprising of DQ/IQ/OQ                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |



## Terms and Conditions:

**Prices** : Prices offered are net-ex-works Mumbai exclusive of Packing, Carting, Crating, Forwarding, Freight and any other incidental charges will be extra to your account.

**Taxes** : **GST (18%) shall be applied.**

**Packing** : **2.5% Extra is Packing And Forwarding Charges**

**Excise** : **NIL**

**Inspection** : Final Inspection prior to dispatch will have to be carried out jointly by us and your representative at our Mumbai works only. You will be intimated as soon as the machines are ready for dispatch.

**Transport** : **Extra to be paid by you directly to the transporter.**

**Delivery** : **Within 16-20 Weeks.**

**Payment** : 1. 40% Advance on confirmation of order.  
2. Balance against proforma invoice through D.D.

**Validity** : This offer is valid for 30 days.

**Warranty** : We undertake to make good by replacement or repair, defects arising out of faulty design, materials or workmanship within 12 months of the date of dispatch, purely on ex-works basis i.e., the defective part will be returned by the buyer to our works and replacement taken by the buyer from our works. Bought-out components are guaranteed by us only to the extent of guarantees given to us by our suppliers. Electrical components, such as heaters, motors, contactors etc. rubber components and instruments such as pressure gauges, thermometers, etc. and other such consumables, etc. are however not covered under this Warranty. Warranty stands void if spare parts other than OEM are used.

**Jurisdiction** : Any dispute arising out of this quotation is subject to the **JURISDICTION OF MUMBAI COURT ONLY**. All dimensions are approximate.

We reserve the right of any technical deviations or modifications in design of all our products due to changes in suppliers situation.

**Thanking you and awaiting for your favorable reply.**

Yours truly,  
**For Hugopharm Technologies Pvt. Ltd.**

Mr. Hiren Panchal  
Tel: 022-24226882 / 022-24222815