



UNICOATER
AUTOMATIC TABLET COATER
USER REQUIREMENT SPECIFICATION

QUOTATION REF NUMBER
SBP/URS/AC/1706/2025-1

JOB NO. : -

REV : 02

CLIENT: **PURPLE LABS PVT LTD**

SHEET NO **1** OF **10**

Equipment: AUTOMATIC TABLET COATER
Document Reference No: SBP/URS/AC/1706/2025-1PURPLE LABS PVT LTD
Contact Person: MR VIJAY CHAUHAN-DIRECTOR

For further assistance please contact:

Sales Office:

Hiren Panchal

S.B. Panchal & Co

8, Jogani Industrial Estate,
541 , Senapati Bapat Marg,
Dadar (W), Mumbai: 400 028
Tel: +91-22-24226882
Fax: +91-22-24222815
Email: sales@sbpanchal.com

REVISION NO.	00					
ISSUED FOR	APPROVAL					
PREPD. BY/DATE	VSP 30/05/25					
CHKD. BY/DATE	VSP 30/05/25					
APPD. BY/DATE	HAP 30/05/25					



**AUTOMATIC TABLET COATER
USER REQUIREMENT SPECIFICATION**

QUOTATION REF NUMBER
SBP/URS/AC/1706/2025-1

JOB NO. : -

REV : 02

CLIENT: **PURPLE LABS PVT LTD**

SHEET NO 2 OF 10

Dear Sir/Madam

Please find below list of options available in standard Automatic Tablet Coater.
Kindly provide the product data in the "Process DataSheet"
Please choose from the various options given in "Vendor Specification Options"
This document shall form the basis of URS and also Design Qualification.

This documents needs to be signed on each page by relevant personnel from Project,
Production, Engineering, Procurement Department.

Sincerely,
For S.B. Panchal & Co.
Hiren Panchal

 HUGOPHARM Systems Engineered with Mind & Spirit	AUTOMATIC TABLET COATER USER REQUIREMENT SPECIFICATION	QUOTATION REF NUMBER SBP/URS/AC/1706/2025-1 JOB NO. : - REV : 02
	CLIENT: PURPLE LABS PVT LTD	SHEET NO 3 OF 10
	INTRODUCTION: Tablet coating systems vary in function and purpose. There are several uses of coaters: • Strengthening coat • De-dust/flavor mask • Color and elegance • Active drug coats • Protective • Timed Release • Sugar coat These systems range from a simple aqueous (water based) coat, to drug layers, and solvent based coating systems. Aqueous systems typically have less complexity, whereas active drug coating layers have a need for more instrumentation and controls, and lastly solvent based systems employ even more instrumentation and controls required primarily for safety reasons (explosion protection). Also, the product itself may pose environmental and exposure hazards that need to be dealt with to implement the installation of the coating system. Lastly, the utilities and space involved needs to be discussed prior to purchasing the equipment. Obviously, the actual equipment size is but a small part of the entire scope of work necessary to implement a larger coating operation. Small coaters can be self-contained, depending upon the application. The unit and its support equipment must be able to be installed in current building facilities. This document encompasses the normal range of coater operations. Technology improvements and new applications may require deviation from this template. This is to be utilized as a guide for the user to answer the majority of questions involved in specifying and using the equipment. Addendums may be used to round out the requirements. A typical Aqueous Process involves the following steps: • Dedust/Preheat (blowers on, pan in jog mode, heat applied, solution recirculating) • After reaching a preset exhaust temperature indicating tablet bed is warm: • Spray Solution to Tablets (Pan in run mode) • After reaching desired weight gain on tablets: • Cooldown (Heat off, Pan in jog mode) • Unload Tablets	

 HUGOPHARM Systems Engineered with Mind & Spirit		AUTOMATIC TABLET COATER USER REQUIREMENT SPECIFICATION		QUOTATION REF NUMBER SBP/URS/AC/1706/2025-1 JOB NO. : - REV : 02
		CLIENT: PURPLE LABS PVT LTD	SHEET NO 4 OF 10	

1. Process Data Sheet: Design Basis		
1.1	Purpose of Coating	• Strengthening coat • De-dust/flavor mask • Color and elegance • Active drug coats • Protective • Timed Release • Sugar coat • Moisture Barrier * Please tick whatever is applicable:
1.2	Activity/Potency of the Product being coated.	The products being used on the machine are of: • Low active (vitamins, supplements) • Medium active (eg. Asprin) • Highly active (eg. Morphine) • Exclusive (eg. Penicillin) * Please tick whatever is applicable:
1.3	Solvent Handled	The coater will be used for (Check all that apply): • Aqueous • Low LEL (Lower Explosion Limit) Solvents • High LEL Solvents • Sugar coat * Please tick whatever is applicable:
1.4	Capacity	The coater shall be capable of producing _____ (unit) of coated tablets per BATCH. OR The coater shall be capable of producing _____ (unit) of coated tablets per 8-hour shift.

 HUGOPHARM Systems Engineered with Mind & Spirit		AUTOMATIC TABLET COATER USER REQUIREMENT SPECIFICATION		QUOTATION REF NUMBER SBP/URS/AC/1706/2025-1
		CLIENT: PURPLE LABS PVT LTD		JOB NO. : - REV : 02 SHEET NO 5 OF 10

Sr. No.	Description	Vendor Specification	User Confirmation
Equipment User Requirement Specification Data Sheet.			
2	Inlet Air Handling Unit		
2.1	Body Type	Single Skinned Double Skinned (Insulated)	
2.2	Material of Construction Inside	S.S. 304 M.S. Painted	
2.3	Material of Construction Outside	S.S. 304 M.S. Painted	
2.4	Primary Filter	20 Micron 10 Micron	
2.5	Secondary Filter	5 Micron 3 Micron	
2.6	Final Filter	0.3Micron HEPA	
2.7	Differential Pressure Across Filters	Gauge Transmitter (On PLC)	
2.8	Heating	Electrical Steam Thermic Fluid	
2.9	Dehumidifier (If product is hygroscopic or effervescent)	Yes No Inlet Air Humidity: • Not critical • Dry at least to _____°C Dew point • Within _____°C to _____°C Dew point	

**AUTOMATIC TABLET COATER
USER REQUIREMENT SPECIFICATION**

QUOTATION REF NUMBER
SBP/URS/AC/1706/2025-1

JOB NO. : -

REV : 02

CLIENT: **PURPLE LABS PVT LTD**

SHEET NO 6 OF 10

2.11	Temperature Control	On/Off	
		PID	
2.12	Air Inlet CFM Indication	Yes	
		No	
2.13	Inlet Air RH Indication	Yes	
		No	
2.14	Air Inlet CFM Control	Yes	
		No	
2.15	Position of Ait Handling Unit	Behind Back Wall	
		On Service Floor	
2.16	Approximate (INSULATED) duct length between AHU and Machine Bottom Chamber		
3	Perforated Pan		
3.1	Material of Construction Inside	S.S. 316	
3.2	Pan Body and Doors Insulation	Yes	
		No	
3.2	Inlet Isolation Valve	On/Off	
		% Opening	
3.3	Interchangeable Bowls Required	Yes (Mention Size)	
		No	
3.4	Bypass Duct Assembly Required for Sugar Coating Intermittent Cycle (Dosing, Rolling, Drying)	Yes	
		No	

AUTOMATIC TABLET COATER USER REQUIREMENT SPECIFICATION

QUOTATION REF NUMBER
SBP/URS/AC/1706/2025-1

JOB NO. : -

REV : 02

CLIENT: **PURPLE LABS PVT LTD**

SHEET NO 7 OF 10

4.	Charging & Discharging	
4.1	Automatic Drum Lift for Charging of Tablets Required	Yes No
4.2	Automatic Unloading of Tablets through Unloading Chute	Yes No
5.	Spraying System	
5.1	Spray Nozzles-AntiBearding Type	Yes No
5.1	Integrated Header required	Yes No
5.2	Spray Pump Motor	Flameproof Non Flameproof
5.3	Solution Tank Heating	Yes (Steam/Electrical) No
5.4	Solution Tank Stirring	Pneumatic Motor Electric Motor
5.5	Spray Rate Indication in PLC (Solution Tank connected to weighing scale)	Yes No
5.6	Drip Nozzles Required for Sugar Coating	Yes No

**AUTOMATIC TABLET COATER
USER REQUIREMENT SPECIFICATION**

QUOTATION REF NUMBER
SBP/URS/AC/1706/2025-1

JOB NO. : -

REV : 02

CLIENT: **PURPLE LABS PVT LTD**

SHEET NO 8 OF 10

6	Exhaust Air System	
6.1	Impeller MOC	Aluminum / M.S. S.S. 304/316
6.2	Blower MOC	M.S. Painted S.S. 304
6.3	Position of Ait Handling Unit	Behind Back Wall On Top Service Floor
6.4	Approximate Duct Length between Machine to Exhaust Blower	
6.7	Exhaust Scrubber	Yes No
7	CIP System	
7.1	CIP Nozzles	Yes No
7.2	CIP High Pressure Pump	Yes No
7.3	Automated Drain Valve in Sink	Yes No
7.4	CIP System with Tank	Yes No

**AUTOMATIC TABLET COATER
USER REQUIREMENT SPECIFICATION**

QUOTATION REF NUMBER
SBP/URS/AC/1706/2025-1

JOB NO. : -

REV : 02

CLIENT: **PURPLE LABS PVT LTD**

SHEET NO 9 OF 10

8	Automation & Control	
8.1	Control Philosophy	Manual PLC SCADA
8.2	PLC Controls	Mistubishi
8.3	Touchscreen Size Preference	5.7" Mitsubishi 10.1" Mitsubishi 12" Mitsubishi
8.4	HMI Based 21CFR Compliance	Yes No
8.5	SCADA and Industrial PC Based 21CFR Compliance.	Yes No
8.6	Print Out	Yes No
8.7	Data Storage on Memory Card/USB on HMI	Yes No Data can only be viewed and printed from HMI
8.8	Data Storage on Remote PC Using Software	Data can only be viewed and printed from the Computer.
8.9	Print Parameters	List of Print Parameter Required
8.10	Flameproof Design for All Motors	Yes No

**AUTOMATIC TABLET COATER
USER REQUIREMENT SPECIFICATION**

QUOTATION REF NUMBER
SBP/URS/AC/1706/2025-1

JOB NO. : -

REV : 02

CLIENT: **PURPLE LABS PVT LTD**

SHEET NO 10 OF 10

10	Certifications-LOCAL COMPLIANCE FOR EXPORTS		
10.1	CE Compliance	Yes No	Cost can increase 30% standard equipment and Lead Time 4 Months
10.2	CE Certified	Yes No	Cost can increase 50% of standard equipment and Lead Time 6 Months
10.3	Atex Compliant	Yes No	Cost can increase 2X standard equipment and Lead Time 6Months
10.4	Atex Certified	Yes No	Cost can increase 3X standard equipment and Lead Time 15-18Months
10.5	Area Classification of Equipment	Zone1 Zone2 Safe Area	